

QUICK SETUP GUIDE

SECULIFE HITAM





The measurement and testing of measuring equipment in the field of medical technology necessitates absolute reliability and a broad range of applications. The SECULIFE HIT_{AM} has been custom tailored to fulfill these demands: It's the ideal device for testing, training and production in the field of medical technology.

The SECULIFE HIT_{AM} TRMS multimeter offers an entirely new dimension of measuring convenience: ultramodern ergonomics which make measuring even easier and more reliable. A broad range of functions that fulfill the requirements of demanding medical technology. And of course world-class quality which is highly convincing all the way down to the finest details.



- 23 multimeter functions
- Voltage measurement
- Auto-ranging current measurement from 100 A (resolution: 10 nA) to 1 A (16 A)
- Capacitance and resistance measurement, diode and continuity testing
- Measuring categories: 600V CAT III and 300V CAT IV
- 1 kHz low-pass filter
- TRMS AC and AC+DC, 20 kHz bandwidth
- Data storage for more than 15,000 measured values
- Housing and protective rubber holster with antimicrobial properties

Interfaces:

- Bidirectional infrared interface for communication with the PC (38.4 kBd)
- Optionally available IR-USB adapter

Power supply: battery or mains operation via optional broad range variable power pack (95 to 250 V AC)

- IP 52 protection against dust and water
- Automatic blocking sockets (ABS)
- Separate battery and fuse compartments
- Made in Germany
- Furnished with DAkkS calibration certificate



Voltage Measurement

- 1. In accordance with the voltage to be measured, turn the rotary switch to V= or V \sim =
- 2. Connect the measurement cables as shown. The connector jack should be grounded. Make sure that a current measuring range ("A") has not been activated when the multimeter is connected for voltage measurement! If the fuse's blowing limits are exceeded as a result of operator error, both the operator and the instrument are in danger! With the rotary switch in the V position, the multimeter is always set to the 1 V measuring range immediately after it is switched on. As soon as the MAN / AUTO key is pressed, and assuming the measured value is less than 90 mV, the multimeter is switched to the mV measuring range.



Alternating Voltage Measurement with 1 MΩLoad Resistance and Frequency Measurement with Selectable Low-Pass Filter

- 1. In accordance with the voltage to be measured, turn the rotary switch to $V_{\sim 1M\Omega}$ or 1kHz
- 2. Connect the measurement cables as shown. The "L" should be gounded

Voltage measurement:

- 1. You can switch back and forth between voltage measurement with and without low-pass filter
- 2. Press the FUNC I ENTER multifunction key repeatedly until the unit of measure V or. V/Fil appears at the display

Frequency Measurement:

- 1. Apply the measured quantity in the same way as for voltage measurement
- 2. Manually select the measuring range for the voltage amplitude. When the instrument is switched to frequency measurement, the previously selected voltage measuring range remains active.
- 3. You can switch back and forth between frequency measurement with and without low-pass filter. Press the **FUNC I ENTER** multifunction key repeatedly until unit of measure Hz or Hz/Fil appears at the display



Connection of the measurement cables

Alternating Voltage and Frequency Measurement V AC and Hz with Selectable Low-Pass Filter

- 1. In accordance with the voltage or frequency to be measured, turn the rotary switch to V~ or Hz.
- 2. Connect the measurement cables as shown. The "⊥" connector jack should be grounded

Voltage measurement:

- 1. You can switch back and forth between voltage measurement with and without low-pass filter
- 2. Press the FUNC I ENTER multifunction key repeatedly until unit of measure V or V/Fil appears at the display

Frequency Measurement:

- 1. Apply the measured quantity is the same way as for voltage measurement
- 2. Manually select the measuring range for the voltage amplitude. When the instrument is switched to frequency measurement, the previously selected voltage measuring range remains active
- 3. You can switch back and forth between frequency measurement with and without low-pass filter. Press the **FUNC I ENTER** multifunction key repeatedly until unit of measure Hz or Hz/Fil appears at the display.



Connection of the measurement cables

Resistance Measurement "Ω"

- 1. Disconnect supply power from the electrical circuit of the device to be measured, and discharge all high-voltage capacitors
- 2. Make sure that the device under test is voltage-free. Interference voltages distort measurement results!
- 3. Set the rotary switch to $_{n}\Omega^{"}$.
- 4. Connect the DUT as shown



Continuity Test

- 1. Disconnect supply power from the electrical circuit of the device to be measured, and discharge all high-voltage capacitors
- 2. Make sure that the device under test is voltage-free. Interference voltages distort measurement results!
- 3. Set the rotary switch to the speaker sign (as shown as in the picture)
- 4. Connect the conductor path under test as shown



Diode Testing with a Constant Current of 1 mA

- 1. Disconnect supply power from the electrical circuit of the device to be measured, and discharge all high-voltage capacitors
- 2. Make sure that the device under test is voltage-free. Interference voltages distort measurement results!
- 3. Set the rotary switch to ->I
- 4. Press the **FUNCI ENTER** key
- 5. Connect the DUT as shown



Temperature Measurement

Measurement with Thermocouples, Temp TC:

- 1. Set the rotary switch to "Temp_{TC}"
- 2. The reference temperature is measured at the internal reference junction
- 3. Connect the sensor to the two accessible jacks. The instrument displays the measured temperature using the selected unit of measure



Measurement with Resistance Thermometers

Set the rotary switch to "Temp_{TC}" or "Temp_{RTD}"

Automatic Compensation:

- 1. Press the **ZERO I ESC** key. "Short leads" appears on the display If you prefer to enter cable resistance directly, you can skip the following entry prompt
- 2. Short circuit the measuring instrument's connector cables. "000.00" appears at the display. After pressing the **FUNC I ENTER** key, automatic compensation of cable resistance is activated for all subsequent measurements. The short-circuit can now be liminated, and the device is ready for use

Entering Cable Resistance:

- 1. Press the **ZERO I ESC** key once again in the automatic compensation menu
- 2. Enter the known resistance of the connector cables with the scroll keys: Select the digit to be changes with the <> keys, and change the respectively selected digit with the up, down keys. The default value is 0.43 Ω . Values can be selected within a range of 0 to 50 Ω
- 3. Upon pressing the **FUNC I ENTER** key, the selected value is activated and the display is returned to the measuring function. Cable resistance remains in memory even after the instrument has been switched off



Capacity Measurement

- Disconnect supply power from the electrical circuit of the device to be measured, and discharge all high-voltage capacitors Make sure that the device under test is voltage-free. Capacitors must always be discharged before measurement is performed. 1.
- 2.
- Set the rotary switch to " \dashv " 3.
- 4. Connect the (discharged!) device under test to the sockets with the measurement cables as shown



Current Measurement

Notes Regarding Current Measurement:

• The multimeter may only be operated with installed batteries or rechargeable batteries. Dangerous currents are otherwise not indicated, and the instrument may be damaged

• Set up the measuring circuit in a mechanically secure fashion, and secure it against inadvertent breaks. Select conductor cross-sections and lay out connections such that they do not overheat

An intermittent acoustic signal warns of current greater than 10 A. A continuous acoustic signal warns of current greater than 16 A

• The input for the current measuring range is equipped with a fuse link. Maximum permissible voltage for the measuring circuit (= rated voltage of the fuse) is 1000 V AC/DC. Use specified fuses only! The fuse must have a breaking capacity of at least 30 kA

• If the fuse for the active current measuring range blows, "FUSE" appears at the digital display, and an acoustic signal is generated at the same time

• If a fuse should blow, eliminate the cause of overload before placing the instrument back into service!

• Be absolutely certain that the measuring ranges are not overloaded beyond their allowable capacities.

Direct and Pulsating Current Measurement, Direct Connection, A DC and A (DC+AC) :

- 1. First disconnect supply power from the measuring circuit or the power consumer (1), and discharge any capacitors
- 2. In accordance with the current to be measured, turn the rotary switch to A= or $A\sim=$
- Select the current type appropriate for the measured quantity by briefly pressing the FUNC I ENTER multifunction key. Each time the key is pressed, the instrument is switched back and forth between A DC and A (DC + AC) TRMS, which is indicated by means of an acoustic signal. The current type is indicated at the LCD by means of the DC or the (DC+AC)TRMS symbol
- 4. Safely connect the measuring instrument (without contact resistance) in series to the power consumer (2) as shown
- 5. Switch supply power to the measuring circuit back on (3)
- 6. Read the display. Make a note of the measured value if the instrument is not being operated in the memory mode or the transmission mode
- 7. Disconnect supply power from the measuring circuit or the power consumer (1) once again, and discharge any capacitors
- 8. Remove the test probes from the measuring point and return the measuring circuit to its normal condition



Alternating Current and Frequency Measurement, Direct Connection, A AC and Hz

- 1. First disconnect supply power from the measuring circuit or the power consumer (1), and discharge any capacitors.
- 2. In accordance with the current or frequency to be measured, turn the rotary switch to A~ or Hz
- 3. Select the desired measured quantity by briefly pressing the **FUNC I ENTER** multifunction key. Each time the key is pressed, AC_{TEMS} and Hz are alternately selected, and switching is acknowledged with an acoustic signal.
- 4. Safely connect the measuring instrument (without contact resistance) in series to the power consumer as shown.
- 5. Switch supply power to the measuring circuit back on (3)
- 6. Read the display. Make a note of the measured value if the instrument is not being operated in the memory mode or the transmission mode
- 7. Disconnect supply power from the measuring circuit or the power consumer (1) once again, and discharge any capacitors
- 8. Remove the test probes from the measuring point and return the measuring circuit to its normal condition



GMC INSTRUMENTS



GMC-I Messtechnik GmbH

Südwestpark 15 \Box 90449 Nürnberg \Box Germany TEL +49 911 8602-111 \Box FAX +49 911 8602-777 www.gossenmetrwawatt.com \Box info@gossenmetrwawatt.com